

WHY YOU CHOOSE US:

- ★ Your satisfaction, value for money and regular references is our motto.
- ★ National & International standard and accreditation for quality and safety is our commitment.
- ★ Long lasting performance & after sales service round the clock is our oath.
- ★ Transparent dealing, ethical business practices is our principle.
- ★ Applying latest technology, regular R&D and innovation to be the market leader & trend setter is our vision.
- ★ Well coordinated team of dedicated professional is our strength.

TECHNICAL SPECIFICATIONS

MODEL	GRAD 50	GRAD 60	GRAD 100
Radiography Rating	50mA, 80 kVp	15-60mA stations 45KVP to 100 KVP in 15 Step	25-100mA stations 45KVP to 100 KVP in 15 Step
Output	3.8 KW	5KW	8KW
Rectification	Full wave rectified	Full wave rectified	Full wave rectified
Timer	0.06 to 5 Sec.	0.06 to 5 Sec.	0.06 to 5 Sec.
Digital Display	-----	Digital Display of mAs, KV & mA	Digital Display of mAs, KV & mA
X-Ray Tube	XD-3/ Imported	XD-3/Imported	XD-3/Imported
Power Supply	230 Voltage AC, 15 Amps.	230 Voltage AC, 15 Amps.	230 Voltage AC, 15 Amps.
Line Regulation	10% 0.6 max	10% 0.4 max	10% 0.4 max
Protection	Electronic over-load protection	Electronic over-load protection	Electronic over-load protection
Tube Stand	Portable / Mobile	Spring balance/ Counter Balanced Mobile	Spring balance/ Counter Balanced/ Mobile/Floor to ceiling Column stand

GRAD MEDICAL EQUIPMENTS PVT. LTD.

Factory : E-93, Sector - 7, Noida (U.P.)
 info@gradmedicalequipments.com
 Website : www.gradmedicalequipments.com
 Mobile : +91 8527760020, Toll Free No. : 1800-313-5700

Marketed by :



Grade in Radiology and Diagnostics

Grad Mobile / Portable Series Line Frequency (50mA, 60mA & 100mA)



GRAD MEDICAL EQUIPMENTS PVT. LTD.



ABOUT US

Grad as its name means grade in radiology & diagnostics is not a saying only but it took 20 years of dedication, hard work, experience & researches to be here as trend setter & upcoming market leader in the radiology field. This is the only company that driven & govern by our esteemed customers & their valuable feedback.

We the Grad is passionate manufacturer specializing in radiology equipments., are manufacturing the excellent quality of products which are designed & developed using most advanced technology & production techniques.

Our consultant & R&D team is working 24x7x365 to understand the requirement, issue, need of client's short term or long term goal & presenting & explaining their findings to the management. These efforts are reflecting in our products & services with a total commitment to quality, safety & customer satisfaction.

PRODUCT RANGE

- Digital Radiographic System
- High Frequency Machines From 4KW to 80 KW Capacity
- Line Frequency Machines from 20mA to 800mA
- Veterinary X-Ray
- Mammography System
- C-Arm Image Intensifier
- Dental OPG Systems
- Radiation Protection Equipments

"Our Motivation is to satisfy exceed customer's demand & expectation and focusing on how to save time, money and nerves in the process of a setup. Its is not enough that we are highly committed and qualified, due to our many years of experience we know and understand what you need and how we can arrange a successful future together."

MOBILE X-RAY UNIT

- Effective design to give excellent and easy operator and patient positioning.
- Digital display of kV, mA, mAs,
- Ruggedly build spring balanced cantilever tube arm
- Easy and safe parking of tube head
- Easy to move in hospital for operation theatre & ward X-Rays.
- Hand Brake system enables easy and effortless movements on slopes, inclines and Ramps
- For hospital room, ICUs, OTs, wards disaster relief, emergency medical treatment, sports events, charity organizations, overseas voyage, etc.



STATE OF ART RESEARCH RESULTS

- High Performance and Reliable Technology
- Less Dose & Best Image Quality
- Sturdy, Low Maintenance & Long Life
- Real Portable Box type Machine
- With dismantle type stand