


YOUR RADIATION SAFETY FOR X RAY FACILITY ROOM WITH

LEAD GLASS

Lead glass is known as x-ray glass or radiation shielding glass which has its calcium component replaced by lead.

GRLGS



LEAD PROTECTION SCREEN (Single Fold/Three Fold)

A lead barrier between x-ray source and operator working on x-ray console protects from radiation. Available single fold 6 feet x 3 feet or in three fold 6 feet x 3 feet joined with 1.5x6 size two additional panel to each side of main panel, to provide full scattered radiation protection. Body mildsteel or wooden ply board with compulsory lead window.



WALL MODEL APRON STAND

The multi chromed swivel caliper having 180° radius of action folded against the wall parallel to each other, the brackets. It saves space.

GRWS



MOBILE SWIVEL ARM APRON STAND

Ergonomic aprons stand for 6 to 12 lead apron on a mobile and can bring anywhere as per need.

GRMS



VIEW BOX

Grad view boxes are made of LED lights behind a Translucent screen and used to provide backlighting for a radiographic image.

GRVB



Grvb may be for single / Double / Triple or four film with automatic film activation sensor and variable brightness controllers.



GRLC

LEAD CURTAIN

Lead Curtain are ideal for areas where secondary or low level radiation is present. We tailored as per our customer requirement.



ABOUT US

Grad as its name means grade in radiology & diagnostics is not a saying only but it took 20 years of dedication, hard work, experience & researches to be here as trend setter upcoming market leader in radiology field. This is the only company that driven & govern by our esteemed customers & their valuable feedback. We the Grad is passionate manufacturer specializing in radiology equipments, are manufacturing the excellent quality of products which are designed & developed using most advanced technology & production techniques. Your satisfaction, value for money and regular references is our motto. National & International standard and accreditation for quality and safety is our commitment. Long lasting performance of our products is our oath. Transparent dealing, ethical business practices is our principle. Applying latest technology, regular R&D and innovation to be the market leader & trend setter is our vision.

GRAD MEDICAL EQUIPMENTS PVT. LTD.

Factory : E-93, Sector - 7, Noida (U.P)
 info@gradmedicalequipments.com
 Website : www.gradmedicalequipments.com
 Mobile : +91 8527760020, Toll Free No. : 1800-313-5700

Marketed by :



Grade in Radiology and Diagnostics

Grad Total Radiation Safety Solution



GRAD MEDICAL EQUIPMENTS PVT. LTD.



RADIATION SAFETY APRONS SERIES

GRAD LEAD APRON FULL BODY WRAP AROUND

GRLWA



These apron are designed to protect your body from back to front through loop and hook



GRLSS

GRAD LEAD APRON SINGLE SIDE

Light Weight Wear Light Easily Wear on and off Comfortable to wear for long durations



SKIRT AND VEST TYPE

GRLSV

Skirt and Vest separation design fully alleviate The pressure off the upper body



GRLFA

LEAD FREE APRONS

We offer lead free aprons, these aprons are environment friendly and very light in weight and crack resistance.



PRODUCT RANGE

- Digital Radiographic System.
- High Frequency Machines From 4KW to 80 KW Capacity.
- Line Frequency Machines from 20mA to 800mA.
- Veterinary X-Ray.
- Mammography System.
- C-Arm Image Intensifier.
- Dental OPG Systems.
- Radiation Protection Equipments.



COMPLETE YOUR RADIATION SAFETY WITH...

LEAD GLOVES :

As hands mostly come active radiation field during radiology procedure so protect your hand from protection.



GRLGV

LEAD CAP :

Lead cap is useful to save occupational head exposure from radiation that may cause of brain cancer.



GRLCP

LEAD GOGGLES :


Grad range of Lead goggle provides the total safety of eyes.



GRLGG

LEAD ACRYLIC FULL FACE SHIELD :

Mask include a ratcheting head strap that gifts all sizes and is especially useful for people who wear glasses and to protect their entire face area from x-ray.



GRFM

THYROID SHIELD :

Due to high Sensitivity to radiation thyroid gland needs special protection. Our GRTS will give you the complete protection to the neck and sternum.



GRTS

GONALD SHIELD MALE :

Gonald Shield for male are the best protection of their reproductive organ.



GRGSM

GONALD SHIELD FEMALE :

Female gonald shield protects their reproductive organs.



GRGSF

LEG PROTECTION SHIELD :

Provides whole protection from radiation with complete legs mobility.



GRLP

BREAST LEAD PROTECTION SHIELD :

Breast lead protection shields most sensitive organs such as Breast, Lungs, Stomach, Colon and bone marrow.



GRBP

FOOT PROTECTION SHIELD :


To Save your foot always wear foot protection shield during radiology procedure.



GRFP

ARMS PROTECTION SHIELD :

Provides arms protection from radiation with complete hand mobility.



GRAP

CUSTOMIZE YOUR ORDER WITH...

SIZE AND LENGHT.

COLOR AND PRINTING.

LED EQUIVALENCE (.25MMPB TO .5MMPB).

STRAP OR VELCRO TYPE.