



Our Motivation is to satisfy exceed customer's demand & expectation and focusing on how to save time, money and nerves in the process of a setup. Its is not enough that we are highly committed and qualified, due to our many years of experience we know and understand what you need and how we can arrange a successful future together.



Grade in Radiology and Diagnostics

## Grad Pink Series mammoGrad (Mamography)

### TECHNICAL SPECIFICATIONS

Model	mammoGrad
X-ray Generator	High Frequency 40kHz
Power output	3.5 KW
Rad kVp Range	20-40kVp
X-ray tube	Rotating anode dual focus type
Rad mAs Range	upto 500 mAs
Operating console	Feather touch digital display of kVp, mA & mAs with APR facility
Interlocks	Voltage error Earthing error kVp error Filament error Thermal error Rotor error
Tube stand	U arm motorized breast compression. Auto & manual shut off. CR compatible. Mammo specialized Cassette size- 18x24 cm & 24x30 cm (as per customer's requirement)
Power requirement	230V, 1 phase AC with line regulation of +/-5%

Due to constant upgradation, company reserves the right to alter the mechanical & technical specs without prior notice.



### GRAD MEDICAL EQUIPMENTS PVT. LTD.

Factory : E-93, Sector - 7, Noida (U.P.)  
 info@gradmedicalequipments.com  
 Website : www.gradmedicalequipments.com  
 Mobile : +91 8527760020, Toll Free No. : 1800-313-5700

Marketed by :




## GRAD MEDICAL EQUIPMENTS PVT. LTD.



## ABOUT US

Grad as its name means grade in radiology & diagnostics is not a saying only but it took 20 years of dedication, hard work, experience & researches to be here as trend setter upcoming market leader in radiology field. This is the only company that driven & govern by our esteemed customers & their valuable feedback. We the Grad is passionate manufacturer specializing in radiology equipments, are manufacturing the excellent quality of products which are designed & developed using most advanced technology & production techniques. Your satisfaction, value for money and regular references is our motto. National & International standard and accreditation for quality and safety is our commitment. Long lasting performance of our products is our oath. Transparent dealing, ethical business practices is our principle. Applying latest technology, regular R&D and innovation to be the market leader & trend setter is our vision.



MamoGrad with the Improvement of Women health awareness becomes a routine examination device in hospital as well as diagnostic Centre.

MamoGrad integrates best image quality, user friendly operation excellent diagnostic accuracy personalized breast care and beautiful design all in one.

MamoGrad is a vital preventative diagnostic Machine to help women and their doctors detect signs of early breast cancer whether experiencing symptoms or not. Early detection is the best prevention. Mammograms can detect abnormal breast tissue approximately two years before it becomes cancerous.

## PRODUCT RANGE

- Digital Radiographic System.
- High Frequency Machines From 4KW to 80 KW Capacity.
- Line Frequency Machines from 20mA to 800mA.
- Veterinary X-Ray.
- Mammography System.
- C-Arm Image Intensifier.
- Dental OPG Systems.
- Radiation Protection Equipments.



MemoGrad is designed with a high level of resistance to the environment factor can be experienced in a mobile facility, This provides the uses with a high level of reliability.

As a vehicle based momoGrad version provides the best solution for the screening and of diagnosis of women's breast